


# Procedure for antimicrobial susceptibility testing of *Haemophilus influenzae* by disc diffusion

Confirm identification of isolate. Sub-culture onto **supplemented chocolate agar** and incubate in a **CO<sub>2</sub>- enhanced atmosphere** (CO<sub>2</sub> incubator or candle-extinction jar).



Prepare inoculum

Prepare a **0.5 McFarland** suspension of the bacteria to be tested in sterile saline

Compare prepared suspension with that of the **0.5 McFarland** standard (control) and adjust turbidity as needed with sterile saline or pure culture until correct density is achieved. Suspension must be used within 15 minutes<sup>1,3</sup>.

Inoculate **Haemophilus Test Medium (HTM) [CLSI guidelines] OR Mueller Hinton agar + 5% defibrinated horse blood and 20mg/L β-NAD [EUCAST]** by streaking the plate with a swab multiple times in different directions to ensure even confluent growth. Allow to dry (Maximum time 10 minutes).

Perform **quality control** of medium as appropriate

Perform **quality control** of antimicrobial discs as appropriate.

CLSI <sup>1</sup>	EUCAST <sup>2</sup>
<b><i>Haemophilus influenzae</i></b>	
<i>Haemophilus influenzae</i> ATCC®49766	
<i>Haemophilus influenzae</i> ATCC®49247	
<i>Escherichia coli</i> ATCC®35218	<i>Staphylococcus aureus</i> ATCC®29213

ATCC® is a registered trademark of the American Type Culture Collection, USA

Place discs on plate with sterile forceps/tweezers/disc dispenser. Do not move the discs once they have touched the agar surface.

Allow discs and media to reach **ambient temperature** before inoculating plates and placing discs

Incubate

35°C±2°C; 5% CO<sub>2</sub>; 16-18 hours




Read **quality control** strain zones of inhibition first. If within limits, read test strain. If quality control fails trouble shoot whole process.

Read zone edges at the point showing no growth against a dark background illuminated with reflected light. Measure zone diameter (mm) using ruler/calipers.

Interpret according to the latest **Clinical and Laboratory Standards Institute (CLSI) / European Committee on Antimicrobial Susceptibility Testing (EUCAST)** guidelines.

Record and report findings.

A rapid β-lactamase test yields clinically relevant information earlier than the results of the antimicrobial susceptibility testing so it should be performed as soon as a *H.influenzae* is identified. **Nitrocefin-based test are the preferred method, chromogenic cephalosporin**, which releases a red compound on hydrolysis by a β-lactamase<sup>3</sup>.



- Clinical and Laboratory Standards Institute (2017); Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Information supplement. CLSI document M100-S27
- European Committee on Antimicrobial Susceptibility Testing; Breakpoint tables for interpretation of MICs and zone diameters. Version 7.1, 2017
- Laboratory methods for the diagnosis of meningitis caused by *Neisseria meningitidis*, *Streptococcus pneumoniae* and *Haemophilus influenzae*; WHO manual; second edition (2011)