

# Procedure for antimicrobial susceptibility testing of *Streptococcus pneumoniae* by disc diffusion

Once identification of isolate is confirmed, sub-culture onto **sheep blood agar** and incubate in a **CO<sub>2</sub>- enhanced atmosphere** (CO<sub>2</sub> incubator or candle-extinction jar) for 18-24 hours<sup>1</sup>.

Prepare inoculum

Prepare a **0.5 McFarland** suspension from sheep blood agar<sup>1,2</sup> or **1 McFarland** suspension from chocolate agar<sup>2</sup> of the bacteria to be tested in sterile saline

Compare prepared suspension of the **0.5 McFarland**<sup>1,2</sup> or **1 McFarland**<sup>2</sup> standard (control) and adjust turbidity as needed with sterile saline or pure culture until correct density is achieved. Suspension must be used within 15 minutes<sup>1,3</sup>.

Inoculate **Mueller Hinton agar (MHA)** with **5% sheep blood plate [CLSI] OR MHA + 5% defibrinated horse blood and 20mg/L β-NAD [EUCAST]** by streaking the plate with a sterile swab multiple times in different directions to ensure even confluent growth. Allow to dry (Maximum 10 minutes)

Perform **quality control** of medium as appropriate

Perform **quality control** of antimicrobial discs as appropriate.

| CLSI <sup>1</sup>               | EUCAST <sup>2</sup>                        |
|---------------------------------|--|
| <i>Streptococcus pneumoniae</i> | <i>Streptococcus pneumoniae</i> ATCC®49619 |

ATCC® is a registered trademark of the American Type Culture Collection, USA

Place discs on plate with sterile forceps/tweezers/disc dispenser. Do not move the discs once they have touched the agar surface.

Allow discs and media to reach **ambient temperature** before inoculating plates and placing discs

Incubate





35°C±2°C; 5% CO<sub>2</sub>; 20- 24 hours

Read zone edges at the point showing no growth against a dark background illuminated with reflected light. Measure zone diameter (mm) using ruler/calipers.

Interpret according to the latest **Clinical and Laboratory Standards Institute (CLSI) / European Committee on Antimicrobial Susceptibility Testing (EUCAST)** guidelines.

Record and report findings.

Read **quality control** strain zones of inhibition first. If within limits, read test strain. If quality control fails trouble shoot whole process.



Use a ruler with a handle if calipers are unavailable

1. Clinical and Laboratory Standards Institute (2017); Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Information supplement. CLSI document M100-S27
2. European Committee on Antimicrobial Susceptibility Testing; Breakpoint tables for interpretation of MICs and zone diameters. Version 7.1,2017
3. Laboratory methods for the diagnosis of meningitis caused by *Neisseria meningitidis*, *Streptococcus pneumoniae* and *Haemophilus influenzae*; WHO manual; second edition (2011)